


SmartVest® SIZING GUIDE

AIRWAY CLEARANCE SYSTEM

This sizing guide will help identify appropriate fitting for the SmartVest® and SmartVest Wrap®. All products are designed to fit a wide range of patients. Each is adjustable to accommodate patients who are slightly smaller or larger than the measurements indicated.

Questions? Call 888-966-2525



CHEST CIRCUMFERENCE

(cm)	(in)	SmartVest®	SmartVest Wrap®
41	16"	Child S-Small	
46	18"	Child S-Small	
51	20"	Child Small	S-Small
56	22"	Child Small	S-Small
61	24"	Child Medium	
66	26"	Child Medium	Small
71	28"	Child Medium	Small
76	30"	Child Large	
81	32"	Child Large	Medium
86	34"	Child Large	Medium
91	36"	Adult Small	
97	38"	Adult Small	Large
102	40"	Adult Medium	Large
107	42"	Adult Medium	Large
112	44"	Adult Large	
117	46"	Adult Large	
122	48"	Adult X-Large	X-Large
127	50"	Adult X-Large	X-Large
132	52"		
137	54"		
152	60"		XX-Large
168	66"		

Note: An Extender may be used on any adult size SmartVest® and will add up to 14 inches to the chest circumference.

The SmartVest Wrap® is intended for patients finding it difficult to fit or apply the standard SmartVest®. Fit is based on the patient's chest circumference at the broadest point.

ELECTROMED, INC.

Creating superior care through innovation®

502 Sixth Avenue N.W.
New Prague, MN 56071
www.Electromed.com